



TwinStrand Duplex Sequencing™

AML-29 Panel (Human), v1.0

Tube Label:

PAN

AML-29 Panel (Human),
v1.0

Part Number: 06-1002-XX

Pipeline ID: aml29-panel-v1.0

Total Panel Size: 58.6 kb

Number of Genes: 29

Number of Target Intervals: 150

Target Type: mutation hotspots or full
coding regions

Genome Reference: hs38DH (Human)

This acute myeloid leukemia (AML) panel is designed to detect ultra low-frequency somatic mutations in 29 genes recurrently mutated in AML patients. 90%–95% of adult AML patients carry a mutation in at least one panel gene. Many of these genes are also recurrently mutated in clonal hematopoiesis of indeterminate potential (CHIP) and myelodysplastic syndrome (MDS).

AML Panel Genes

<i>ASXL1</i>	<i>FAM5C (BRINP3)</i>	<i>IDH2</i>	<i>NPM1</i>	<i>RAD21</i>	<i>TET2</i>
<i>CBL</i>	<i>FLT3</i>	<i>KIT</i>	<i>NRAS</i>	<i>RUNX1</i>	<i>TP53</i>
<i>CEBPA</i>	<i>GATA2</i>	<i>KRAS</i>	<i>PHF6</i>	<i>SMC1A</i>	<i>U2AF1</i>
<i>DNMT3A</i>	<i>HNRNPK</i>	<i>MLL-X (KMT2A-X)</i>	<i>PTEN</i>	<i>SMC3</i>	<i>WT1</i>
<i>EZH2</i>	<i>IDH1</i>	<i>MYH11-CBFB</i>	<i>PTPN11</i>	<i>STAG2</i>	

Full Target List

Gene	Accession #	Target Region
<i>ASXL1</i>	NM_015338	a.a. 363-1542
<i>CBL</i>	NM_005188	a.a. 366-477
<i>CEBPA</i>	NM_001287424	full coding region
<i>DNMT3A</i>	NM_022552	a.a. 286-913
<i>EZH2</i>	NM_004456	a.a. 87-208, 244-302, 503-752
<i>FAM5C (BRINP3)</i>	NM_199051	a.a. 80-142, 396-767
<i>FLT3</i>	NM_004119	a.a. 569-647, 807-847
<i>GATA2</i>	NM_032638	a.a. 77-481
<i>HNRNPK</i>	NM_002140	a.a. 21-85, 173-215, 319-336, 371-453
<i>IDH1</i>	NM_005896	a.a. 106-138
<i>IDH2</i>	NM_002168	a.a. 126-178
<i>KIT</i>	NM_000222	a.a. 412-448, 788-828
<i>KRAS</i>	NM_004985	a.a. 1-96



Full Target List - Continued		
Gene	Accession #	Target Region
<i>MLL-X (KMT2A-X)</i>	NM_005933 (<i>KMT2A</i>)	<i>MLL</i> intron 9
<i>MYH11-CBFB</i>	NM_022844 (<i>MYH11</i>)	<i>MYH11</i> intron 30, exon 31
<i>NPM1</i>	NM_002520	a.a. 258-282, 283-295
<i>NRAS</i>	NM_002524	a.a. 1-96
<i>PHF6</i>	NM_032458	full coding region
<i>PTEN</i>	NM_001304717	a.a. 258-337, 385-440
<i>PTPN11</i>	NM_002834	a.a. 47-110, 484-533
<i>RAD21</i>	NM_006265	full coding region
<i>RUNX1</i>	NM_001754	full coding region
<i>SMC1A</i>	NM_006306	a.a. 38-99, 447-515, 578-637, 687-732, 772-902, 1096-1145
<i>SMC3</i>	NM_005445	a.a. 184-268, 365-435, 656-705, 882-1035, 1100-1158
<i>STAG2</i>	NM_001282418	a.a. 42-128, 155-297, 436-472, 513-546, 578-675, 787-844, 892-1155
<i>TET2</i>	NM_001127208	a.a. 1-2003
<i>TP53</i>	NM_000546	full coding region
<i>U2AF1</i>	NM_006758	a.a. 16-44
<i>WT1</i>	NM_024426	a.a. 372-523

a.a. = amino acids

Panel-Specific Protocol Recommendations

Cot-1 DNA Recommended:	Human: 03-2020-XX	Number of Final PCR Cycles:	14 (PCR-2), 6 (PCR-3)
Number of Captures Required:	1	Conditioning Reagent:	Recommended
Typical On-Target Percentage:	80%-90% (1 capture), 90%-95% (2 captures)		

Panel-Specific Protocol Notes

DNA Conditioning: The Library Conditioning step is recommended for this panel. Refer to the TwinStrand Duplex Sequencing Kit Manual for a full explanation of the benefits vs. costs of performing Library Conditioning.

1 vs. 2 Captures: With panels of this size, it is possible to perform only a single round of hybrid capture. However, some AML-29 panel targets have pseudogenes or nearby repetitive elements that can affect efficiency of capture and sequencing. A single round of hybrid capture is generally sufficient to obtain 80% or greater on-target bases with this panel. We recommend starting with 2 captures the first time through the protocol to generate data that can be compared against. An increase in off-target reads means that raw sequencing must be increased proportionally to generate the same amount of informative data, so a trade-off exists between time saved with a single capture vs. potential additional sequencing cost required. The user should consider the pros and cons in the context of their specific situation.

Final PCR Cycles: The number of final PCR cycles varies with probe panel, number of captures, input DNA mass, input DNA quality, thermocycler calibration and other factors. If 2 rounds of hybrid capture are performed then PCR-3 is the final PCR, and with 1 round of capture PCR-2 is the final PCR. The number of final PCR cycles listed in the table for 1 vs. 2 captures should be considered a starting point, assuming 500 ng of high quality genomic DNA input at the beginning of library preparation.

Statistical Power: To detect ultra-low frequency mutations, sufficient Duplex Sequencing Depth is required. Duplex Depth is the denominator when calculating mutation frequency. For example, to observe a mutation frequency of 1/5,000 a library must be sequenced to at least 5,000x Duplex Depth across the specific genomic locus of interest. Additionally, since DNA molecules are being randomly sampled from a pool, to have 95% power to detect a mutation at a frequency of 1/X, a library should be sequenced to at least 3X depth at the locus of interest. For example, to achieve 95% power for detecting a mutation with a true frequency of 1/5,000, a depth of 15,000x is required at that locus. Although TwinStrand™ probe panels are designed for maximum capture and sequencing depth uniformity, coverage cannot ever be fully uniform due to natural sequence context variation. Some sites of interest may achieve below average depths and this should be accounted for in power calculations.



DNA Input vs. Projected Duplex Sequencing Data					
DNA Input (ng)*	Number of Clusters (Million)**	Number of Paired-End Reads (Million)**	Illumina Sequencing Cost (NovaSeq, S4, List Price)***	Mean Duplex Depth****	Max Duplex Depth****
250	75	150	\$216	5,500x	8,000x
500	150	300	\$434	11,000x	16,000x
1,000	300	600	\$864	22,000x	32,000x

* Assumes high molecular weight genomic DNA (quantified pre-fragmentation). DNA requirements to achieve a given Duplex Depth when using cell-free DNA or heavily damaged DNA such as that from FFPE will be lower and higher, respectively.

** Approximate sequencing required to achieve a peak tag family size ~10.

*** S4 list price = \$7,200/lane (2.5 billion clusters/lane) as of December 2019.


**** Data estimates assume 150 bp paired-end reads

i TwinStrand Duplex Sequencing data yields are dependent on accurate quantification, and the quality of input DNA. We recommend Qiagen silica column-based kits or “salting out” methods for DNA extraction. Phenol-chloroform extraction can damage DNA and reduce data yields. DNA should be quantified using a double-stranded DNA-binding fluorescent dye (e.g. Qubit fluorometer, Thermo Fisher Scientific) rather than UV spectrometry which can overestimate concentration if contaminants are present, which will reduce data yields. If possible, check genomic DNA quality using a 0.5% agarose gel, or a fragment analyzer system (e.g. Agilent TapeStation system). DNA integrity number (DIN) >7 on a fragment analyzer or a single high molecular weight band on a gel generally indicates good quality DNA, although oxidation or other forms of non-fragmentation chemical damage may still be present and not appreciable.

Notes



Notes



©2020 TwinStrand Biosciences, Inc. All rights reserved. All trademarks are the property of TwinStrand Biosciences, Inc. or their respective owners. This Information Sheet may contain trademarks, service marks and trade names of other parties, which are the property of their respective owners. Our use or display of third parties' trademarks, service marks, trade names or products herein is not intended to, and does not imply a relationship with, or endorsement or sponsorship by them or us.

